

Vaccine Storage Temperature Log (Fahrenheit)

Location of Unit: _____ Month/Year: _____

- **Twice a day** place an "X" in the box that corresponds with the current temperature and time. If the temperature is not on the grid, write the actual temperature in the "OTHER" box; do not use an X.
- Fill in the minimum and maximum temperatures at the same time (if you have that information).
- Take immediate action on any temperature you find in the shaded area!
- If the temperature in your refrigerator/freezer is not on the grid **or** is in a shaded area, mark the vaccine in question and do not use it until you've called the Minnesota Immunization Program, 800-657-3970 or 651-201-5522.
- If the vaccine has gotten too warm or too cold, you will not be able to distinguish viable from nonviable vaccine by looking at it.
- Keep these monthly temperature logs for 3 years.

Refrigerator - Aim for 40°F

Days 1-16

	DAY	1		2		3		4		5		6		7		8		9		10		11		12		13		14		15		16		
		TIME	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm		
Too warm	min/max																																	
	OTHER																																	
	49°	Take immediate action if temperature is in shaded area Adjust thermostat and call MDH for guidance																																
	48°																																	
47°																																		
Store at 35° - 46°F. Aim for 40°F!	46°																																	
	45°																																	
	44°																																	
	43°																																	
	42°																																	
	41°																																	
	40°																																	
	39°																																	
	38°																																	
	37°																																	
	36°																																	
35°																																		
Too cold	34°	Take immediate action if temperature is in shaded area Adjust thermostat and call MDH for guidance																																
	33°																																	
	OTHER																																	

Freezer - Aim for 0°F

	DAY	1		2		3		4		5		6		7		8		9		10		11		12		13		14		15		16	
		TIME	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	
Too warm	min/max																																
	OTHER																																
	8°	Take immediate action if temperature is in shaded area Adjust thermostat and call MDH for guidance																															
	7°																																
6°																																	
Store at 5°F or ↓	5°																																
	4°																																
	3°																																
	OTHER																																
	-58° or ↓																																

