

# Supply Checklist for Point Prevalence Survey (PPS)

## COVID-19 PPS TOOLKIT FOR LONG-TERM CARE FACILITIES

Consult with testing laboratory to determine capacity and requirements for collection materials, storage, forms, and shipping. The following will be needed for specimen collection.

- ☐ Test kits (swabs)
- ☐ Specimen collection tubes with saline or viral transport media
- ☐ Specimen tube labels, pre-written with resident information
- ☐ Biohazard zip lock bags to store/transport specimens (provided by laboratory)
- ☐ Laboratory forms, completed before testing day
- ☐ Paper, pencil, and thin sharpies for labelling sample tubes (if needed)
- ☐ Transport container/box, ice packs (enough to keep all specimens cool as directed by laboratory). Specimens should be stored at 2–8°C (35–47°F)
- ☐ Trash bags
- ☐ Cart to transport materials through facility, one cart per Testing Team
- ☐ Printed or electronic lists of residents and staff and pens, as needed
- ☐ Symptom assessment form (if desired)

---

## Testing Team Supplies

All staff should practice universal masking with medical-grade face masks and universal front and side eye protection (e.g., face shield, goggles).

The following is needed for each team member collecting specimens assisting.

- N95 respirators or medical-grade face masks (if respirators are not available).
  - One respirator or face mask for each person involved in specimen collection, plus extras in case of contamination or malfunction.

## SUPPLY CHECKLIST FOR POINT PREVALENCE SURVEY (PPS)

- Face shields.
  - One face shield for each person involved in specimen collection. If soiled during specimen collection, should be disinfected with wipes.
- Gowns.
  - If supplies allow, disposable gowns can be changed between every person being tested. If supplies are limited, gowns can be reused if no resident/patient care is being provided. Gowns should always be disposed of if soiled (e.g., visible debris, cough or sneezed upon) or if they come into direct contact with a resident.
- Gloves.
  - Must be changed between every person being tested.
  - Ensure availability of one box, appropriately sized, for each Swabber.
  - Swabbing Assistants also need appropriately sized gloves.
- Alcohol-based hand rub.
- Disinfectant wipes.



Minnesota Department of Health | [health.mn.gov](http://health.mn.gov) | 651-201-5000  
625 Robert Street North PO Box 64975, St. Paul, MN 55164-0975

Contact [health.communications@state.mn.us](mailto:health.communications@state.mn.us) to request an alternate format.

06/24/2020